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



## Handling procedure for Sintered Discs

**Heraeus**

**EUROPE**  
**Heraeus Quarzglas**  
**GmbH & Co. KG**  
 Base Materials  
 Reinhard-Heraeus-Ring 29  
 D-63801 Kleinostheim  
 Phone +49 (6181) 35-7444  
 Fax +49 (6181) 35-7210  
 E-Mail sales.basematerials@heraeus-quarzglas.com

**USA**  
**Heraeus Quartz America LLC**  
 HQI Base Materials Division  
 2400 Grand Avenue Parkway  
 Suite 102  
 Austin, TX 78728  
 Phone +1 (512) 989-05 03  
 Toll free 866 - 989-05 05  
 Fax +1 (512) 989-05 04  
 E-Mail sales@heraeus-quartz.com

**ASIA**  
**Shin-Etsu Quartz Prod. Co., Ltd.**  
 Shinjuku San-Ei Bldg., 12th Floor  
 22-2, 1-Chome, Nishi Shinjuku  
 Shinjuku-Ku  
 Tokyo 160  
 Phone +81 (333) 48 19 11  
 Fax +81 (333) 48 49 19  
 E-Mail d3sales@sqp.co.jp

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## Handling procedure for Sintered Discs

### Temperature Change, Drying and Sterilizing

Sudden temperature change and uneven heating should be avoided. Sintered glass filter funnels, pipeline filter tubes and other sintered glass apparatus with disc diameters exceeding 50 mm, which are to be dried or sterilized, should be placed in cold ovens or sterilizers. The rate of heating should not be more than 2 °C/min. This is the only way of preventing internal strains which are caused by excessive temperature differences between the surrounding glass vessel and the sintered filter disc, and which can lead to fracture of apparatus.

Filtration apparatus should, whenever possible, stand on its rim (stem upwards) in the oven or sterilizer. A perforated support base is advantageous since it allows air convection between the inside of the vessel and the body of the oven. If angled positioning of the filtration apparatus in the oven is unavoidable (pipeline filter tubes), then any point of support which is near to the filter seal position should be protected against premature heatings. This is done by using an underlay of heat-insulating material.

The apparatus should remain in the oven or sterilizer during cooling. Due to the thermal inertia of this type of oven, the cooling time is adequate.

### Cleaning new Sintered Glassware

Before using sintered glass filter apparatus for the first time, hot hydrochloric acid followed by several rinses of distilled water should be sucked through the filter disc under a good vacuum. This removes dust particles and powdered glass. It is important that each successive water rinse be started only after the preceding one has been completely flushed through. This so-called "tear through" method must only be used for cleaning filters. It should never be adopted for preparative or analytical filtration.

<http://www.heraeus-quarzglas.com>

Printed in Germany - Technical data and pictures given in this brochure are subject to change

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# Handling procedure for Sintered Discs

## Mechanical Cleaning

*Sintered glass filters should always be cleaned immediately after use*

If no precipitate has entered the pores, surface rinsing under the tap or with a wash bottle is often sufficient. The filter disc surface can be wiped clean with a small brush or squeegee.

Where some of the precipitate has entered the pores, back-flushing is necessary. In the case of porosities 0 to 2 this can be done simply by using a water tap, connecting it with rubber tubing to the stem of the sintered glass piece and allowing water to run backwards through the filter disc. The water pressure must not exceed 1 kp/cm<sup>2</sup>.

For porosities 3, 4 and 5 the precipitate is flushed or wiped off the disc, and water is sucked through in the opposite direction to filtration.

Filters clogged by dust and dirt during gas filtration can be restored by treatment with a warm detergent solution followed by blowing through clean air from the clean side of the filter.

Dirt particles are brought to the surface by the foam and removed by rinsing with water.

## Chemical Cleaning

If, after mechanical cleaning, some of the pores still remain clogged, or if it is desirable to make sure that no residue from previous work remains before filtering a new substance, then thorough chemical cleaning is required. The choice of solvent obviously depends on the nature of the contamination. For example:

### Barium sulphate

hot conc. sulphuric acid (100 °C)

### Silver chloride

hot ammonia liquor

### Red copper oxide

hot hydrochloric acid and potassium chlorate

### Mercury residue

hot, conc. nitric acid

### Mercury sulphide

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### Albumen

hot ammonia liquor or hydrochloric acid

### Grease, oil

carbon tetrachloride

### Other organic substances

hot, conc. sulphuric acid with an addition of nitric acid, sodium nitrate, or potassium dichromate

### Animal charcoal

careful heating with a mixture of 5 vol. of conc. sulphuric acid + 1 vol. conc. nitric acid to about 200 °C

Prolonged rinsing with water must obviously follow

For biochemical work, cleaning with dichromate sulphuric acid should be avoided, since trivalent chromium compounds, present or newly formed by reduction, are absorbed on the surface of the filter disc. When they are released during subsequent use, biological substances can be seriously damaged. This danger is eliminated by using sulphuric acid with a nitrate or perchlorate addition. Only easily soluble reduction products are formed which can be completely removed by re-washing with water. Since hot, concentrated phosphoric acid and hot alkaline solutions attack the glass surface, they are unsuitable as cleaning agents. If they have to be filtered, an increase in pore size and thus reduced life of the apparatus is unavoidable.


Fig. 1: Water and air flow rate through filter discs of various porosities as a function of pressure differential. For disc diameter 30 mm.